



**Subject: Mask's Color Clarification**

The purpose of this letter is to clarify the subject of the Mask Color. As stated in the Mask specifications, portion of mask's the raw material is made of Copper powder. This raw material has an impact on the Mask color; hence, color variations can be expected in different manufacturing lot.

I.e., the raw copper material used in different production batches, may affect the mask's color, and can vary in the range of orange-brown-copper color (Pic 1) to Salomon pink color (pic 2).

We would like to state that the mask Antiviral and Antibacterial properties are not affected by the mask's color.

**Pic 1 – Orange-Brownish**



**Pic 2 – Salmon-Pink**



Sincerely

**Alon Raivit**

A handwritten signature in black ink, appearing to read 'Alon Raivit', written over a white background.

CEO

[www.olvitech.co.il](http://www.olvitech.co.il) | [info@olvitech.co.il](mailto:info@olvitech.co.il)